IFU


(Instruction for use)

Análogo / Analogue / Análogo

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The analog is part of the family of laboratory components and is used as a replica of the installed implant in the patient, to perform laboratorial work.

This cervical part can simulate several types of fittings: morse taper, external hexagon, internal hexagon and abutments. In its apical part, it has a geometry for axial and radial retention to avoid any type of movement in the working model. The raw material used is surgical stainless steel 300 series, according to ASTM F 899 303.

The analogues are presented in the table below:

Sales Code	Connection Type							Diameter	Connection Measure	Rosca
	CM	LZ	HE	Transmucoso						
				Multi	Mini 4,0-5,0	Mini 4,8	Angul.			
9780275			X					3,5	2,7	M2
9780274			X					4,1	2,7	M2
9780276			X					5,0	2,7	M2
9780538	X	X						3,5	Hex. 2,0	M1.6
9780539	X	X						4,0	Hex. 2,0	M1.6
9780540	X	X						5,0	Hex. 2,0	M1.6
9780540						X	X	4,8	Rotac.	M1.4
9780343				X				4,0	Hex. 2,95	M2
9780344				X				5,0	Hex. 2,95	M2

APPLICATION

Analogues are indicated for simulating the dental implant connection and position in a working model in the laboratory phase. Before handling in the laboratory, make sure that the elements are properly attached.

CONTRAINDICATIONS

We do not recommend reuse, this procedure may directly bad influences to the quality and success of the procedure.

HANDLING

After obtaining the transfer template, screw the analog to the respective impression transfer. For screw tightening, use manual torque with the bi-digital key. Pour the plaster over the mold covering all the analog. After drying the plaster, remove the transfer template from it, thus obtaining the study model or the working model.

PRESENTATION

Any analogue is indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied ready to use because they undergo a decontamination process that eliminates any kind of dirtiness.

STERILIZATION

Supplied non-sterile. Indicated for single and laboratorial use.

PRECAUTIONS

Do not use if the product packaging is broken;

Check the correspondence between the diameter of the analog and the diameter of the installed implant;

Beware of handling, keeping the integrity of the platform, avoiding damage and foreign bodies at the docking interface and screw hole;

After the prosthesis is made, evaluate the appropriate adaptation over the analogue;

Installation torques higher than recommended or use of improper connection keys, can cause damage to the material and make the system unusable;

Do not reuse analogue. That procedure can cause damage to the interface, compromise the positioning and the screw fixing;

Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty;

The component must be perfectly seated to the analog before applying the necessary torque on the screw.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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