IFU



(Instruction for use)

Base de Titânio Ti-Base 4 CEREC® / Titanium Base Ti-Base 4 CEREC® / Base de Titanio Ti-Base 4 CEREC®

PORTUGUÊS

ENGLISH

ESPAÑOL

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CÓDIGO DE VENDAS	DESCRIÇÃO	INTERFACE PROTÉTICA		Ø IMPLANTE (MM)	ALTURA DE CINTA (MM)	PARAFUSO DE FIXAÇÃO	
		CM	HE			ROSCA	CÓDIGO
9780612	TI-BASE 4 CEREC® CM 0,8MM	X		3,5-4,0-5,0	0,8	M1,6	9780545
9780613	TI-BASE 4 CEREC® CM 1,5MM	X		3,5-4,0-5,0	1,5	M1,6	9780545
9780614	TI-BASE 4 CEREC® CM 3,0MM	X		3,5-4,0-5,0	3,0	M1,6	9780545
9780615	TI-BASE 4 CEREC® CM 4,5MM	X		3,5-4,0-5,0	4,5	M1,6	9780545
9780501	TI-BASE 4 CEREC® HE Ø3,5		X	3,5	0,6	M2,0	9780281
9780502	TI-BASE 4 CEREC® HE Ø3,5 X 1,5MM		X	3,5	1,5	M2,0	9780281
9780503	TI-BASE 4 CEREC® HE Ø3,5 X 3,0MM		X	3,5	3,0	M2,0	9780281
9780504	TI-BASE 4 CEREC® HE Ø3,5 X 4,5MM		X	3,5	4,5	M2,0	9780281
9780505	TI-BASE 4 CEREC® HE Ø4,1		X	4,1-5,0	0,6	M2,0	9780281
9780506	TI-BASE 4 CEREC® HE Ø4,1 X 1,5MM		X	4,1-5,0	1,5	M2,0	9780281
9780507	TI-BASE 4 CEREC® HE Ø4,1 X 3,0MM		X	4,1-5,0	3,0	M2,0	9780281
9780508	TI-BASE 4 CEREC® HE Ø4,1 X 4,5MM		X	4,1-5,0	4,5	M2,0	9780281
9780616	SCAN POST 4 CEREC® CM	X		3,5-4,0-5,0	5,2	M1,6	9780545
9780509	SCAN POST 4 CEREC® HE		X	3,5-4,1-5,0	5,2	M2,0	9780281

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This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results, it should be used only with original associated components According to AllPrime Select instructions and recommendations. All procedures must be performed by qualified professionals with specific training in CAD/CAM technology, in the Sirona Dental CAD/CAM System for digitizing and in the production of custom prostheses.

DESCRIPTION

TITANIUM BASE TIBASE 4 CEREC® AllPrime Select: The Titanium Base TiBase 4 CEREC® AllPrime Select is an intermediate prosthetic abutment Indicated for unitary cementation or bolted restorations on dental implants installed in the jaw. It is manufactured in Titanium alloy, according to ASTM F136, to be installed between the implant and the prosthesis, which is manufactured from the Sirona Dental CAD/CAM System for digitizing and production of custom prostheses. The Titanium Base TiBase 4 CEREC® AllPrime Select Has a cylindrical shape with indexator guide (Notch) to fit the ScanBody, prosthesis or Coping. It has a cemable area height of 4.7 mm and with variations in height of Transgingival (see table). The geometry of the cemable area of the titanium Base TiBase 4 CEREC® AllPrime Select is compatible with the blocks Incoris, the Zi Meso Or Disilicate Available for CEREC®. The Titanium Base TiBase 4 CEREC® AllPrime Select is available for Use In implants AllPrime Select of Interfaces Morse Taper and External Hexagon. The fixing screws do not accompany the titanium bases TiBase 4 CEREC® AllPrime Select, are sold separately. Note: Single use Only.

SCANPOST 4 CEREC® AllPrime Select: Scanpost is an extender device about implants For scanning, with plug equivalent to the TiBase 4 CEREC® AllPrime Select that can be used Intra-oral Or Extraorally (Benchtop scanner). It is manufactured in titanium alloy, according to ASTM F136 Standard, has a cylindrical shape with indexator guide (Notch) to fit the ScanBody And is bolted In Implant to the prosthetic connection. The ScanPost 4 CEREC® AllPrime Select can be used up to 50 times. If the screws are lost or damaged, they must be replaced. The ScanPost 4 CEREC® AllPrime Select is available for use in AllPrime Select implants of Morse Taper and External Hexagon Interfaces. The fixing screws do not accompany the ScanPost 4 CEREC® AllPrime Select, are sold separately.

⚠ Note: the Scanpost 4 CEREC® AllPrime Select should not be used as a definitive prosthetic abutment!

⚠ Note: For greater scanning accuracy, AllPrime Select recommends the acquisition of Original Scanbodies Dentsply Sirona.

SALES CODE	DESCRIPTION	PROSTHETIC INTERFACE		Ø IM-PLANT (MM)	TRANSMU-COSAL HEIGHT (MM)	SCREW	
		MT	EH			THREAD	CODE
9780612	TI-BASE 4 CEREC® CM 0,8MM	X		3,5-4,0-5,0	0,8	M1,6	9780545
9780613	TI-BASE 4 CEREC® CM 1,5MM	X		3,5-4,0-5,0	1,5	M1,6	9780545
9780614	TI-BASE 4 CEREC® CM 3,0MM	X		3,5-4,0-5,0	3,0	M1,6	9780545
9780615	TI-BASE 4 CEREC® CM 4,5MM	X		3,5-4,0-5,0	4,5	M1,6	9780545
9780501	TI-BASE 4 CEREC® HE Ø3,5		X	3,5	0,6	M2,0	9780281
9780502	TI-BASE 4 CEREC® HE Ø3,5 X 1,5MM		X	3,5	1,5	M2,0	9780281
9780503	TI-BASE 4 CEREC® HE Ø3,5 X 3,0MM		X	3,5	3,0	M2,0	9780281
9780504	TI-BASE 4 CEREC® HE Ø3,5 X 4,5MM		X	3,5	4,5	M2,0	9780281
9780505	TI-BASE 4 CEREC® HE Ø4,1		X	4,1-5,0	0,6	M2,0	9780281
9780506	TI-BASE 4 CEREC® HE Ø4,1 X 1,5MM		X	4,1-5,0	1,5	M2,0	9780281
9780507	TI-BASE 4 CEREC® HE Ø4,1 X 3,0MM		X	4,1-5,0	3,0	M2,0	9780281
9780508	TI-BASE 4 CEREC® HE Ø4,1 X 4,5MM		X	4,1-5,0	4,5	M2,0	9780281
9780616	SCAN POST 4 CEREC® CM	X		3,5-4,0-5,0	5,2	M1,6	9780545
9780509	SCAN POST 4 CEREC® HE		X	3,5-4,1-5,0	5,2	M2,0	9780281

The abutments 4 CEREC® AllPrime Select are presented in the table below:

⚠ Note: For External Hexagon implants diameter 5,0mm, the use of platform is indicated Switching.

APPLICATION

TITANIUM BASE TIBASE 4 CEREC® AllPrime Select: The

Titanium Base TiBase 4 CEREC® AllPrime Select is indicated according to the space Interocclusal Available, height of transmucosal and three-dimensional position of the implant. They should be used in cementation-based unitary restorations on the titanium Base TiBase 4 CEREC® AllPrime Select which are screwed on the implant. All interfaces are indicated for bolted unit restorations on dental implants. Check the compatibility between the chosen prosthetic interfaces. Titanium Bases are available in the Morse Taper (MT) and External Hexagon (EH) prosthetic interfaces, and are compatible with implants with the respective prosthetic interfaces. It is recommended to use the technique CAD/CAM technology through the Sirona Dental CAD/CAM System.

SCANPOST 4 CEREC® AllPrime Select: The Scanpost 4 CEREC® AllPrime Select is indicated to facilitate scanning in hard to reach areas. The ScanPosts are available in the Morse Taper (MT) and External Hexagon (EH) prosthetic interfaces, and are compatible with implants with the respective prosthetic interfaces. It is recommended to use the technique CAD/CAM technology through the Sirona Dental CAD/CAM System.

CONTRAINDICATIONS

- The Titanium Base TiBase 4 CEREC® AllPrime Select does not allow customization;
- It is contraindicated in restorations with angulation correction of more than 20 ° in relation to the implant axis and restorations whose length exceeds a ratio of 1:1.25 to Implant length (Example: Implant 10mm, not exceed 12,5mm);
- It is contraindicated for cases of multiple restorations;
- Not compatible with the Software Sirona Galileos®
- This product is contraindicated for patients presenting signs of allergy or hypersensitivity to the chemical components of the material: Titanium alloy Ti6Al4V-ELI (ASTM F136);

MAGNETIC RESONANCE IMAGING (MRI) – INFORMATION

These products are manufactured in metallic material that can be affected by MRI energy.

For more information, see “Security Information | Magnetic resonance imaging (MRI) “in <https://www.allprimedental.com>

HANDLING

FOR INTRAORAL SCANNING:

TI BASE 4 CEREC® AllPrime Select: Select and install the Titanium Base TiBase 4 CEREC® AllPrime Select About the Implant according to previous planning. Enter the Transfer Scanning (Scanbody)* Compatible With the equipment of the CAD/CAM system used on the titanium Base TiBase 4 CEREC® AllPrime Select, and then run the Scanning according to manufacturer’s instructions Sirona. Make sure the correct fitting of the Scanbody* About the Titanium Base TiBase 4 CEREC® AllPrime Select.

⚠ Note: After selecting the NBB 3.4 L Library, in the module of choice, select Tibase.

SCANPOST 4 CEREC® AllPrime Select: Select and install The Scanpost 4 CEREC® AllPrime Select About the Implant according to previous planning. Enter the Transfer Scanning (Scanbody)* Compatible With the CAD/CAM system equipment used on the The Scanpost 4 CEREC® AllPrime Select, and then run the Scanning according to manufacturer’s instructions Sirona. Make sure the correct fitting of the Scanbody* About the Scanpost 4 CEREC® AllPrime Select.

⚠ Note: After selecting the NBB 3.4 L Library, in the module of choice, select ScanPost.

* the Transfer Scanning (Scanbody) compatible with titanium Base Tibase 4 CEREC® AllPrime Select is an accessory supplied by Dentsply Sirona.

FOR GYPSUM MODEL SCANNING:

TI BASE 4 CEREC® AllPrime Select: Perform the molding of the implant according to the Previous planning, patient’s clinical situation and implant interface. Use artificial gums In the plaster model to simulate the soft tissue of the patient. Manufacture the plaster model according With the appropriate techniques. Select and install the Titanium Base TiBase 4 CEREC® AllPrime Select On the implant analog. Enter The Transfer Scanning (Scanbody)* Compatible with the equipment of the CAD/CAM system used and perform the scanning according to the manufacturer’s instructions. Make sure the correct Fitting the Transfer On the titanium Base TiBase 4 CEREC® AllPrime Select.

⚠ Note: After selecting the NBB 3.4 L Library, in the module of choice, select TiBase.

SCAN POST 4 CEREC® AllPrime Select: Perform the molding of the implant according to the Previous planning, patient’s clinical situation and implant interface. Use artificial gums In the plaster model to simulate the soft tissue of the patient. Manufacture the plaster model according With the appropriate techniques. Select and install Scanpost 4 CEREC® AllPrime Select on the implant analog. Enter the Transfer Scanning (Scanbody)* Compatible with the equipment of the CAD/CAM system used and perform the scanning as shown in the of the manufacturer. Make sure the correctFitting the Transfer About the Scanpost 4 CEREC® AllPrime Select.

⚠ Note: After selecting the NBB 3.4 L Library, in the module of choice, select ScanPost.

* the Transfer Scanning (Scanbody) compatible with titanium Base Tibase method 4 CEREC® AllPrime Select is an accessoriesupplied by Dentsply Sirona.

MESOSTRUCTURE DESIGN:

Design the prosthesis structure in Sirona InLab Software (version 3.65) or Sirona CEREC® Software (version 4.2) According to previous prosthetic planning and the material to be used. Titanium Bases TiBase 4 CEREC® AllPrime Select About implants AllPrime Select (MT and EH) Can be used in conjunction with the products Sirona, presented in the following table:

Note: the ScanPost Compatibility according to the table below:

⚠ Note: The Software does not permit the use of Scan Posts 4 CEREC® AllPrime Select for jobs where Ti-Bases 4 CEREC®

SIRONA PRODUCTS FOR USE IN CONJUNCTION WITH TIBASE 4 CEREC® AllPrime Select				
LIBRARY	SCANBODY (SIRONA)	REF. SCANBODY OMNICAM (SIRONA)	REF. SCANBODY BLUECAM inEos X5 and inEos Blue (SIRONA)	CEREC BLOCKS RECOMMENDED
NBB 3.4L	L	6431329	6431303	inCoris ZI meso F0.5 (L) inCoris ZI meso F2 (L) IPS e.max CAD (L)

AllPrime Select are used with 1,5mm / 3,0mm e 4,5mm transmucosal heights.

PRODUCTS FOR USE IN CONJUNCTION WITH SCANPOST 4 CEREC® AllPrime Select					
SCANPOST 4 CEREC® (AllPrime Select)	LIBRARY	SCANBODY (SIRONA)	REF. SCANBODY OMNICAM (SIRONA)	REF. SCANBODY BLUECAM inEos X5 and inEos Blue (SIRONA)	TIBASE 4 CEREC® (Me-dens)
9780616	NBB 3.4L	L	6431329	6431303	9780612
9780509					9780501
					9780505

⚠ Note: The thickness of the minimum bone wall around the wall of the screw grooving of the structure is variable according to the material, as shown in the following table. The maximum angulation of the structure shall not exceed 20°. The concavity of the structure shall not exceed 6° (when Coping). In the case of an angled structure the height of the cembale portion from the emergency profile (prosthetic height), should not exceed 1:1.25 in ratio to Implant length.

MATERIAL	MINIMUM THICKNESS(MM)
ZIRCONIA	0,5
IPS e-max® CAD	0,9

- Be sure not to exceed a 20 ° angle between the implant shaft and the shaft of the restoration.
- If the Mesostructural is designed to receive the aesthetic ceramics, make sure that this will not narrow the screw grooving.
- The fitting cavity in the titanium Base Tibase 4 CEREC® AllPrime Select mesostructural should not be coated.
- Make sure there are no sharp edges or corners in the design of the Mesostructural.
- Before making the structure of the prosthesis, make sure that the size of the machining block is compatible with the design to be carried out.
- Manufacture the workpiece as shown, using machining equipment Sirona CEREC® MC X, InLab MC XL or InLab MC X5.

CLINICAL PROOF OF PROSTHESIS STRUCTURE:

1. Perform a clinical proof of the structure Made on the titanium Base Tibase 4 CEREC® AllPrime Select.
2. Make sure the proper fitting between The parts clinically and through radiography.
3. Make sure the occlusal space exists, Proximal Inter And that the structure of the prosthesis meets Aesthetic and functional requirements.
4. Finish the restoration according to the planning Appropriate techniques.

CEMENTING OF THE TITANIUM BASE TIBASE 4 CEREC® AllPrime Select:

1. It is recommended to screw the titanium Base Tibase 4 CEREC® AllPrime Select On an implant analog during The entire finish and finishing process of the prosthesis.
2. Position the Titanium Base Tibase 4 CEREC® AllPrime Select On the implant analog and transcend the Fastening on the same with a slight tightening Using the key Installation indicated (see table below).
3. The surfaces of the titanium Base Tibase 4 CEREC® AllPrime Select Aimed at cementing the structure of the Prosthesis should be treated with sandblasting (Aluminum oxide, 50 µm, maximum pressure of 2 bar). Protect screw access during The cementing process.
4. ⚠ Note: It is recommends the use of chemically resinous cements Activated for bonding in metals (e.g.: PANAVIA – Kuraray). In Disilicate Lithium (IPS E.max® CAD) is mandatory the use of IVOCCLAR cement Multilink. The handling of the cement must follow the Manufacturer’s instructions.
5. Apply Cement to External portion of the titanium Base Tibase 4 CEREC® AllPrime Select and press the reset byfitting according With the Indexer tab.
6. The restoration must be Pressed on the titanium Base Tibase 4 CEREC® AllPrime Select And the excess of cement that is extravasated by the Orifice of the same shouldbe immediately Removed.
7. After the cement’s prey Infrastructure is unscrewed

from implant analog and Remove excess cement remaining in the Titanium Base Tibase 4 CEREC® AllPrime Select.

CUSTOM PROSTHESIS INSTALLATION:

Sanitizer and sterilize the set, titanium Base Tibase 4 CEREC® AllPrime Select and custom prosthesis, before To takeit in the mouth. After these procedures, position the setin the mouth using the installation driver and the torque indicated as shown in the following table. Perform the occlusal adjustments and Interproximal needed.

SCREW MODEL	TORQUE (N.cm)	INSTALLATION DRIVER
9780545	25	HEXAGONAL 1,2MM
9780281	30	HEXAGONAL 1,2MM

PRESENTATION FORM

All parts that make up the system 4 CEREC® AllPrime Select are individually packaged in surgical paper packaging type Tyvek® and transparent film display. They are provided clean because they undergo a decontamination process, which eliminates any kind of dirt, however, require sterilization prior to clinical use.

STERILIZATION

The Tibase 4 CEREC® AllPrime Select are supplied Non-sterile. All models are indicated for single use only.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization by Saturated Steam - Autoclaves

In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.

All surgical instruments and prosthetic components must be sterilized prior to use.

Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if The packaging is violated;
- Check the correspondence between The type of screw according to the implant model with its prosthetic interfaces;
- Before installing the products, make sure that they have the appropriate form of use for the selected prosthetic system;
- Special care must be taken during the project of the prosthesis so that the space Interocclusal No Insufficient or excessive.
- Ensure that the parts are not swallowed or aspirated by the person who is receiving the product (patient);
- All models of TiBases 4 CEREC® AllPrime Select Are Indicated For single use only. Reuse of this product can Cause:

- *Effects Biological adverse effects resulting from residual residues, be they, microorganisms and/or substances arising from previous uses;*

IFU | Instruction for use

- *Changes in physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may impair the functionality of which it was intended;*
- *Reuse This product does not guarantee its safety and efficacy, and exempts any warranty from the related products.*

- Higher insertion torques than recommended or use of improper installation driver can cause damage to the material and make the system unusable;
- Always use the product sequence AllPrime Select. The use of instrumental and/or Third-party components do not Ensures the perfect function of the sYstem of iMplante AllPrime Select and exempts any warranty from the product.
- Improper planning can compromise the performance of the set Prosthesis, resulting in system failures such as loss or fracture of the implant or component, looting or fracture of the prosthetic screw;
- At the time of installation Of the component in the implant, Transgingival or implant analog, be sure to align The connection and/or o Screw, thus avoiding damage to the thread or locking of the system;
- The Titanium Base Tibase 4 CEREC® AllPrime Select is not indicated for cementation in the mouth;
- Communicate to the laboratory that any modification or damage to the prosthetic interface requires that the workpiece be discarded and remade;
- When choosing titanium Base Tibase 4 CEREC® AllPrime Select Ensure that it is compatible with the interface and diameter of the implant installed;
- The Titanium Base Tibase 4 CEREC® AllPrime Select is not suitable for multiple restorations;
- Check the compatibility of Scanbodies, the With the Sirona software used;
- The components of the system 4 CEREC® AllPrime Select Must Be perfectly seated to implant After applying the required torque to the Screw. This check is recommended using Periapical radiographs;

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

CEREC® is a registered trademark of Sirona Dental Systems GmbH.

Galileos is a mark used by the Sirona Dental Systems GmbH Germany.

IPS E.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.

Multilink® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.

Panavia™ is a registered trademark of Kuraray Co, LTD, JP.

The aforementioned libraries are protected by copyright and belong to their respective authors. All rights reserved to such authors. These libraries are available in Sirona InLab Software (version 3.65) or Sirona CEREC® Software (version 4.2) of the Sirona Dental Systems GmbH. All rights reserved.

Not all products are available in all countries. Favor Contact Your authorized Distributor.

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CÓDIGO DE VENTAS	DESCRIPCIÓN	INTERFAZ PROTÉSICA		ØIMPLANTE (MM)	ALTURA DEL TRANSMUCOSO (MM)	TORNILLO DE SUJECIÓN	
		MT	EH			TORNILLO	CÓDIGO
9780612	TI-BASE 4 CEREC® CM 0,8MM	X		3,5-4,0-5,0	0,8	M1,6	9780545
9780613	TI-BASE 4 CEREC® CM 1,5MM	X		3,5-4,0-5,0	1,5	M1,6	9780545
9780614	TI-BASE 4 CEREC® CM 3,0MM	X		3,5-4,0-5,0	3,0	M1,6	9780545
9780615	TI-BASE 4 CEREC® CM 4,5MM	X		3,5-4,0-5,0	4,5	M1,6	9780545
9780501	TI-BASE 4 CEREC® HE Ø3,5		X	3,5	0,6	M2,0	9780281
9780502	TI-BASE 4 CEREC® HE Ø3,5 X 1,5MM		X	3,5	1,5	M2,0	9780281
9780503	TI-BASE 4 CEREC® HE Ø3,5 X 3,0MM		X	3,5	3,0	M2,0	9780281
9780504	TI-BASE 4 CEREC® HE Ø3,5 X 4,5MM		X	3,5	4,5	M2,0	9780281
9780505	TI-BASE 4 CEREC® HE Ø4,1		X	4,1-5,0	0,6	M2,0	9780281
9780506	TI-BASE 4 CEREC® HE Ø4,1 X 1,5MM		X	4,1-5,0	1,5	M2,0	9780281
9780507	TI-BASE 4 CEREC® HE Ø4,1 X 3,0MM		X	4,1-5,0	3,0	M2,0	9780281
9780508	TI-BASE 4 CEREC® HE Ø4,1 X 4,5MM		X	4,1-5,0	4,5	M2,0	9780281
9780616	SCAN POST 4 CEREC® CM	X		3,5-4,0-5,0	5,2	M1,6	9780545
9780509	SCAN POST 4 CEREC® HE		X	3,5-4,1-5,0	5,2	M2,0	9780281

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CEREC® es una marca registrada de Sirona Sistemas dentales GmbH.
Galileos es una marca utilizada por el Sirona Sistemas dentales GmbH Alemania.
IPS E.max® es una marca registrada de Ivoclar Vivadent AG, Liechtenstein.
Multilink® es una marca registrada de Ivoclar Vivadent AG, Liechtenstein.
Panavia™ es una marca registrada de Kuraray Co, LTD, JP.