IFU

(Instruction for use)

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The healing cups are part of the prosthetic components family of transient use, and have 2 distinct parts: head and body. The head has an internal hexagon of 1.2 mm for coupling in the driver that will allow the torque application, besides different interfaces, according to the type of implant where the component will be installed. The head has different heights that correspond to the gingival heights, taking into account the variations of mucosal thickness and also presents with different diameters or sizes, what is called a regular or esthetic design. The body has a thread with geometry that allows the fitting to the specific implant.

The raw material used is the titanium alloy (grade 5), according to ASTM F136.

The healing cups are presented in the table below:

| SalesCode | Interface Type | | | Implant Diameter | Size | | Height | Head Diameter | Thread | Maximum Torque |
|-----------|----------------|----|----|------------------|---------|----------|--------|---------------|--------|----------------|
| | CM | LZ | HE | | Regular | Esthetic | | | | |
| 9780250 | | | X | 3,5 | X | | 1,5 | 4,0 | M2 | 20 |
| 9780251 | | | X | 3,5 | X | | 3,0 | 4,0 | M2 | 20 |
| 9780252 | | | X | 3,5 | X | | 4,5 | 4,0 | M2 | 20 |
| 9780253 | | | X | 4,0 | X | | 1,5 | 4,1 | M2 | 20 |
| 9780254 | | | X | 4,0 | X | | 3,0 | 4,1 | M2 | 20 |
| 9780255 | | | X | 4,0 | X | | 4,5 | 4,1 | M2 | 20 |
| 9780256 | | | X | 5,0 | X | | 1,5 | 5,0 | M2 | 20 |
| 9780257 | | | X | 5,0 | X | | 3,0 | 5,0 | M2 | 20 |
| 9780258 | | | X | 5,0 | X | | 4,5 | 5,0 | M2 | 20 |
| 9780617 | X | | | 3,5 | X | | 1,5 | 3,5 | M1.6 | 15 |
| 9780618 | X | | | 3,5 | X | | 3,0 | 3,5 | M1.6 | 15 |
| 9780619 | X | | | 3,5 | X | | 4,5 | 3,5 | M1.6 | 15 |
| 9780620 | X | | | 4,0 | X | | 1,5 | 4,0 | M1.6 | 15 |
| 9780621 | X | | | 4,0 | X | | 3,0 | 4,0 | M1.6 | 15 |
| 9780622 | X | | | 4,0 | X | | 4,5 | 4,0 | M1.6 | 15 |
| 9780623 | X | | | 5,0 | X | | 1,5 | 5,0 | M1.6 | 15 |
| 9780624 | X | | | 5,0 | X | | 3,0 | 5,0 | M1.6 | 15 |
| 9780625 | X | | | 5,0 | X | | 4,5 | 5,0 | M1.6 | 15 |
| 9780528 | X | | | 3,5-4,0-5,0 | X | | 1,5 | 3,5 | M1.6 | 10 |
| 9780529 | X | | | 3,5-4,0-5,0 | X | | 3,0 | 3,5 | M1.6 | 10 |
| 9780530 | X | | | 3,5-4,0-5,0 | X | | 4,5 | 3,5 | M1.6 | 10 |
| 9780531 | X | | | 3,5-4,0-5,0 | X | | 6,0 | 3,5 | M1.6 | 10 |

| SalesCode | Interface Type | | | Implant Diameter | Size | | Height | Head Diameter | Thread | Maximum Torque |
|-----------|----------------|----|----|------------------|---------|----------|--------|---------------|--------|----------------|
| | CM | LZ | HE | | Regular | Esthetic | | | | |
| 9780259 | | | X | 3,5 | | X | 1,5 | 4,5 | M2 | 20 |
| 9780260 | | | X | 3,5 | | X | 3,0 | 4,5 | M2 | 20 |
| 9780260 | | | X | 3,5 | | X | 4,5 | 4,5 | M2 | 20 |
| 9780262 | | | X | 4,0 | | X | 1,5 | 4,6 | M2 | 20 |
| 9780263 | | | X | 4,0 | | X | 3,0 | 4,6 | M2 | 20 |
| 9780263 | | | X | 4,0 | | X | 4,5 | 4,6 | M2 | 20 |
| 9780265 | | | X | 5,0 | | X | 1,5 | 5,5 | M2 | 20 |
| 9780266 | | | X | 5,0 | | X | 3,0 | 5,5 | M2 | 20 |
| 9780267 | | | X | 5,0 | | X | 4,5 | 5,5 | M2 | 20 |
| 9780626 | X | | | 3,5 | | X | 1,5 | 4,0 | M1.6 | 15 |
| 9780627 | X | | | 3,5 | | X | 3,0 | 4,0 | M1.6 | 15 |
| 9780628 | X | | | 3,5 | | X | 4,5 | 4,0 | M1.6 | 15 |
| 9780629 | X | | | 4,0 | | X | 1,5 | 4,5 | M1.6 | 15 |
| 9780630 | X | | | 4,0 | | X | 3,0 | 4,5 | M1.6 | 15 |
| 9780631 | X | | | 4,0 | | X | 4,5 | 4,5 | M1.6 | 15 |
| 9780632 | X | | | 5,0 | | X | 1,5 | 5,5 | M1.6 | 15 |
| 9780633 | X | | | 5,0 | | X | 3,0 | 5,5 | M1.6 | 15 |
| 9780634 | X | | | 5,0 | | X | 4,5 | 5,5 | M1.6 | 15 |
| 9780532 | X | | | 3,5-4,0-5,0 | | X | 1,5 | 4,5 | M1.6 | 10 |
| 9780533 | X | | | 3,5-4,0-5,0 | | X | 3,0 | 4,5 | M1.6 | 10 |
| 9780534 | X | | | 3,5-4,0-5,0 | | X | 4,5 | 4,5 | M1.6 | 10 |
| 9780535 | X | | | 3,5-4,0-5,0 | | X | 6,0 | 4,5 | M1.6 | 10 |

APPLICATION

The healing cups are indicated for the creation of the gingival tunnel that will give access to the implant in the prosthetic phase. It can be used on the implant installation or the re-opening surgery (second surgical phase). It is indicated according to the existing gingival height and three-dimensional position of the implant.

Check the compatibility between the fitting of the healing cup and the implant and, at the time of the installation of the healing cup, make sure that the component has been tightened with the necessary torque according to the table above.

CONTRAINDICATIONS

We do not recommend reuse, this procedure may directly bad influences to the quality and success of the procedure.

MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information | Magnetic Resonance Imaging (MR)” at <https://www.allprimedental.com>

HANDLING

After opening the package, deposit the healing cup on a sterile surface and collect it with the 1.2 mm hexagonal driver and install it on the implant of the same interface.

At the time of installation, be sure to align the healing cup line along the implant axis.

The component must be perfectly seated to the implant after the necessary torque has been applied to the healing cup. The maximum torque value to be applied to each wound should be respected. It is recommended that the healing cup is 1mm above the gingival level.

PRESENTATION

Any healing cup is indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization By Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this

temperature is reached.

- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Maximum care is recommended when using the instruments so that they are not swallowed or aspirated by the patient;
- Any type of healing cup is for single use only. Reusing this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- It is recommended that the healing cup is remove from the patient’s mouth within 6 months after installation;
- This product should not be used as a support for removable prostheses;
- Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty;
- The component must be perfectly seated to the implant after the necessary torque has been applied to the healing cup. This check is recommended using periapical radiographs;
- Instruct the patient that if he or she notices that the healing cup will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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