IFU


(Instruction for use)

Compactador e Expansor / Compactor and Expander / Compactador e Expansor

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The AllPrime Select bone compactors form a set of 4 pieces that have a stepped cylindrical shape and have different diameters. At its active tip, the product has laser markings that refer to the implant lengths existing in the AllPrime Select system, so the user can detect the depth of compaction performed for subsequent implant installation. At the end of the active tip, the product has a sharpened concave geometry. At the other tip, the bone compactor has a cylindrical shape with a special fit to be attached to the AllPrime Select hand-driver or even used as an impact tool for the AllPrime Select surgical hammer.

The raw material used in all of the drivers is tempered stainless steel, according to ASTM F899.

The bone compactors models are presented in the table below:

Bone Compactors			
Sales Code	Diameter		Active Length
	Initial	Final	
9780816	2,0	3,0	46
9780817	2,8	3,7	46
9780818	3,2	4,2	46
9780819	4,0	5,0	46

The AllPrime Select bone expanders form a set of 6 pieces that have a cylindrical shape with different diameters. At its active tip, the product has laser markings that reference the existing implant lengths in the AllPrime Select system, so the user can detect the depth of the expansion performed for the subsequent installation of the implant. At the end of the active tip, the product has a convex geometry. At the other tip, the bone expander has a cylindrical shape with a special fit to be attached to the AllPrime Select hand-driver or even used as an impact tool for the AllPrime Select surgical hammer.

The raw material used in all of the drivers is tempered stainless steel, according to ASTM F899.

The bone expanders models are presented in the table below:

Bone Expanders		
Sales Code	Diameter	Active Length
9780820	2,0	46
9780821	2,5	46
9780822	3,0	46
9780823	3,5	46
9780824	4,0	46
9780825	4,5	46

The AllPrime Select threaded bone expanders form a set of 6 pieces that have a cylindrical shape with different diameters. At its active tip, the product contains a thread that facilitates and guides the insertion of the part into the bone. It also has laser markings that reference to the implant lengths existing in the AllPrime Select system, so the user can detect the depth of the expansion performed for the posterior installation of the implant. At the end of the active tip, the product has a convex geometry. At the other tip the threaded bone expander has a universal square fit of 3.9 mm, which connects to the gear of the AllPrime Select torque wrench.

The raw material used in all of the drivers is tempered stainless steel, according to ASTM F899.

The threaded bone expanders models are presented in the table below:

Threaded Bone Expanders		
	Diameter	Active Length
9780809		
9780810	2,0	25
9780811	2,5	25
9780812	3,0	25
9780813	3,5	25
9780814	4,0	25
9780809	4,5	25

APPLICATION

The AllPrime Select bone compactor is a manual surgical instrument used exclusively in the maxilla for the elevation of the maxillary sinus floor, preparing the surgical alveolus for the installation of implants in the posterior maxillary region with little bone height.

The AllPrime Select expander and threaded expander are manual surgical instruments used to perform the bone expansion technique that prepares the surgical alveolus for the installation of implants in cases of low thickness ridge. They are only indicated where the thickness of the alveolar crest has a minimum of 3 mm. The threaded expander can be used as a bone tap in high density mandibular bones.

CONTRAINDICATIONS

The compactors and expanders do not present contraindications, provided they are used correctly for the purposes indicated.

HANDLING

Once sterilized, the compactor or expander should only be handled in a sterile environment by properly trained professionals and in appropriate costume at the time of surgery for the installation of dental implants.

All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the instruments in the best possible way, respecting the AllPrime Select indications for each instrument.

Bone compactor:

Select the bone compactor that meets the implant diameter to be installed. It should be used in the posterior maxilla when there is insufficient bone height to implant installation. With the help of the surgical hammer and the hand-driver, insert the compactor into the

bone cavity carefully, without breaking the sinus mucosa, to the desired length.

Bone expander:

It should be used when the bone thickness is thin, with a minimum of 3 mm.

Threaded expander bone:

It should be used with the help of the universal torque wrench when the bone thickness is thinner with at least 3 mm and requires greater strength.

Note:

- 1 - Always keep the same direction of the instrument during insertion.
- 2 - The pressure on the compactor or expander must be controlled.

SANITATION AND STERILIZATION

The AllPrime Select compactors and expanders are reusable, provided non-sterile and should be sanitized and sterilized before and after use by following the steps mentioned in the link below:

PRESENTATION

The AllPrime Select compactors and expanders are individually packaged in tyvek® surgical paper and transparent film display. The instruments are also sold in kits, assembled in identified position cases, and fixed by silicone rings. The kit is guaranteed by exclusive seals in both the length and width of the case. The instrument undergoes a decontamination process that eliminates any type of dirtiness, but the sanitation process, described above, must be followed as soon as its first used.

PRECAUTIONS

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- During the use of the compactor, be careful not to break the sinus mucosa of the maxillary sinus;
- During the use of the expander, should be careful not to cause fractures in the lateral bone walls of the surgical alveolus;
- Should be careful in the compactor sanitation process because it is a sharp metal instrument, avoiding any type of physical trauma;
- Be sure to use the diameter of the bone compactor or expander that is compatible with the procedure and the diameter of the planned implant;

- The instruments should not be re-sharped;
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be difficult;
- Maximum care is recommended when using the instruments so that they are not swallowed or aspirated by the patient;
- Make sure you have all the necessary instruments to perform the surgery according to surgical planning;
- Before each procedure, check the conditions of the AllPrime Select instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear.
- The use of excessive force can compromise the surgical act;
- Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty.

ADVERSE EFFECTS

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can cause damage to the surgical sequence of the implant to be installed.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

LIFE CYCLE

AllPrime Select compactors and expander should be disposed when it loses its functionality or the user detects excessive wear.

SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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