**IFU**

*(Instruction for use)*

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

## DESCRIPTION

The bone chisels form a set of 6 pieces between short and long models of different diameters, depending on the model of cover screw implant overlapped by bone to be removed. The bone chisel, at the end of its active tip, has a geometry with 4 sharpened tips and a circular guide that fits within the hexagonal socket of the implant cover screw. At the other tip, the bone chisel has an anatomically designed body for bi-digital use, and has an independently-driven washer and hole for dental floss installation, in order to increase safety during use. Although the washer has circular movements independent of the body, it is fixed and can't be disassembled from the instrument.

The raw material used is tempered stainless steel, according to ASTM F899.

The bone chisels are presented in the table below:

Código de Venda	Tamanho		
	Diâmetro de Corte	Comprimento	
		Curto	Longo
9780832	3,7	X	
9780833	4,4	X	
9780834	5,4	X	
9780835	3,7		X
9780836	4,4		X
9780837	5,4		X

## APPLICATION

The bone chisel is a surgical instrument used to remove bone that eventually overlaps the exposed implant cover screw, preventing its bi-digital removal for posterior prosthetic component installation. Usually that can occur with infra-bony installed implants.

## CONTRAINDICATIONS

Bone Chisels do not present contraindications, provided they are used correctly for the purposes indicated.

## HANDLING

Once sterilized, the bone chisel should be used in compliance with all biosafety and surgical asepsis requirements.

All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the instruments in the best possible way, respecting the AllPrime Select indications for each instrument.

Select the bone chisel that meets the installed implant cover diameter, also choose the chisel that has the best length for the procedure, depending on the position of the implant in the mouth.

For patient safety, thread and tie a dental floss into the hole in the washer at the top of the chisel, giving greater safety to the procedure. Take the bone chisel to the point where the implant is installed and already exposed, and perform circular movements over the implant cover until the overlapping bone is removed, allowing the bi-digital release of the implant cover screw.

## SANITATION AND STERILIZATION

The bone chisels are reusable, provided non-sterile and should be sanitized and sterilized before and after use by following the steps

mentioned in the link below:

[Surgical Kit Cleaning Manual](#)

## PRESENTATION

The bone chisels are individually packaged in tyvek® surgical paper and transparent film display. The instrument undergoes a decontamination process that eliminates any type of dirtiness, but the sanitation process described above must be followed as soon as it's first used.

## PRECAUTIONS

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- Should be careful in the sanitation process with the scalpel because it is a sharp metal instrument, avoiding any type of physical trauma;
- Be sure to use the diameter of the bone chisel compatible with the installed implant diameter;
- Bone chisel should not be re-sharpened;
- The non-replacement of a bone chisel, as recommended by AllPrime Select, may not properly perform the removal of the bone overlying the implant cover, compromising the success of the procedure;
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be difficult;

- Maximum care is recommended when using the instruments so that they are not swallowed or aspirated by the patient. Always use the dental floss technique on the washer hole at the top of the chisel, giving greater safety to the procedure;
- Make sure you have all the necessary instruments to perform the surgery according to surgical planning;
- Before each procedure, check the conditions of the AllPrime Select instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear.
- The use of excessive force may compromise the procedure;
- Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty.

## ADVERSE EFFECTS

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can cause damage to the surgical sequence of the implant to be installed.

## STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

## DISPOSAL OF MATERIAL

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

## LIFE CYCLE

This product is recommended for up to 30 uses, provided that AllPrime Select-recommended conditions of use are followed. Regardless of the number of times the instrument has been used, the professional must always evaluate the instrument's conditions after each use.

## SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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