## IFU


(Instruction for use)

# Implante Cone Morse / Morse Taper Implant / Implante Cono Morse

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllSelect Prime instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

**DESCRIPTION**

The AllSelect Prime morse taper implants are made of titanium (grade 5), according to ASTM F136.

It has an external surface mechanically treated with a blasting operation with aluminum oxide, finished with acid attack. In its external macrogeometry the morse taper implant has an apical taper to facilitate the beginning of its insertion in the surgical alveolus, the external thread has variable pitch, which provides high primary stability, besides the double profile of this thread, which confers a higher insertion velocity of the implant. The conical fitting has an angle of 10 °, with a 2.0 mm hexagonal prosthetic indexer, which is also used for installation, besides the M1.6 thread for the components fixing. The morse taper implant has a single prosthetic interface, independent of its diameter, for components exclusively of morse taper connection, but in its cervical region the morse taper implant has a conical seat, exclusively used in components of the zero line, in other words, the AllSelect Prime has two prosthetic lines for the same implant, giving greater versatility to the morse taper implant. The specific implant cover screw has a maximum installation torque of 10N. cm and is provided in the implant package and is located attached to the top of the transparent tube cap. The maximum installation torque for morse taper implants is 70 N.cm.

The AllSelect Prime morse taper implants are presented in the table below:

Morse Taper Implant								
Sales Code	Diameter	Lenght	Sales Code	Diameter	Lenght	Sales Code	Diameter	Lenght
9780515	3,5	8,5	9780521	4,0	8,5	9780527	5,0	8,5
9780510	3,5	10,0	9780516	4,0	10,0	9780522	5,0	10,0
9780511	3,5	11,5	9780517	4,0	11,5	9780523	5,0	11,5
9780512	3,5	13,0	9780518	4,0	13,0	9780524	5,0	13,0
9780513	3,5	14,5	9780519	4,0	14,5	9780525	5,0	14,5
9780514	3,5	16,0	9780520	4,0	16,0	9780526	5,0	16,0

**APPLICATION**

The AllSelect Prime morse taper implant is indicated for surgical procedures in maxilla or mandible bones providing support for prosthetic components. It can be used in one or two-stage processes, for single or multiple prostheses, in an immediate or conventional loading protocol, depending on the primary stability, implant diameter and length, occlusal load and all other biomechanical analyzes required for this decision. In properly planned situations, the implant can be installed immediately after extraction of dental root. According to the literature, implant success is based on obtaining a direct and rigid connection between the living bone and the surface of the implant, from surgery performed under controlled conditions.

Morse taper implants are compatible with the prosthetic components that present the morse taper or zero line prosthetic interface of the AllSelect Prime system. They are not compatible with components of other systems or prosthetic interfaces.

In case of doubts about the indication, the use must be interrupted until all the questionable points have been clarified.

**CONTRAINDICATIONS**

It is not recommended to install the implant in the presence of acute infectious or inflammatory processes, inadequate bone volume or quality, serious medical problems such as disorders of bone metabolism, disorders of blood coagulation, inadequate healing, inadequate oral hygiene, incomplete jaw growth, patient uncooperative and unmotivated, drug or alcohol abuse, psychoses, prolonged functional disorders that resist any drug treatment, dry mouth, weakened immune system, diseases that require the regular use of steroids, endocrinological diseases, pregnancy or any other disorder, be it systemic or psychological, which may harm the tissue repair process.

**MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION**

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information

| Magnetic Resonance Imaging (MR)” at [www.allprimedental.com](http://www.allprimedental.com)

**HANDLING**

The dental implant should be used according

all the requirements of biosafety and surgical asepsis.

The implant should be handled only in a sterile environment, by properly trained professionals and in appropriate attire at the time of surgery. The same must happen with the other instruments that will be used in the procedure.

All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the components and instruments in the best possible way, respecting the AllSelect Prime indications for each product.

**Drilling:**

Once the surgical bed has been made, the implant must be removed from its packaging through the corresponding driver, taken to the artificial alveolus and inserted by threading until the planned length.

Drill sequence recommended by AllSelect Prime in

Morse Taper Implant						
Drill Diameter	Bone Type					
	I and II			III and IV		
	Implant Diameter			Implant Diameter		
	3,5	4,0	5,0	3,5	4,0	5,0
2	X	X	X	X	X	X
2,5	X	X	X	X	X	X
2,8	X	X	X	X	X	X
3,2		X	X	X	X	X
3,5		X	X		X	X
3,8			X		X	X
4,5			X			X
4,8						X

drilling for implant installation:

In the bone bed drilling procedure, it is recommended to check the instruction for use of the AllSelect Prime surgical drills for more details on the application, handling and precautions of the instrument.

**Implant installation:**

Identify the planned implant model to be installed through the label on the back of the carton pack. The identification is given easily through text and colors, in the case of morse taper implants, the color is dark blue. This label also identifies the diameter, length, batch number, date of manufacture and the validity of the implant. The diameter is also easily identified by the colors of the plastic tube cap on the front of the package: white for 3.5 mm diameter implants, red for 4.0 mm diameter implants and black for 5.0 mm diameter implants. At the top of this cap there is a circular label that identifies the length of the implant. In this front part of the package, we can also identify the circular label indicative of sterilization, which will appear in red color, if this label is in another color, or the product

is not sterilized or has the sterilization expired and, in this case, the implant should not be used.

To remove the implant from the packing, remove the carton to access the internal blister and pull the tyvek surgical paper through the corner of the blister, accessing the inner tube containing the implant without the need of gloves. Lay the sterile tube over the surgical field without touching it.

With sterile surgical gloves, hold the bottle with a non-dominant hand and open the lid with the other hand, always keeping the tube upright, that is, with its opening facing up, preventing the implant from leaving the packaging.

In the opening of the cap, the user has access to the implant connection, which can be collected with the specific installation driver, searching the perfect fit between the driver and the implant by pressing lightly against each other. The installation driver has a mechanical collection system that secures the implant to the key, and helps in its removal from the transparent tube of the package, without requiring any other manipulation than your own driver. Transport the implant to the surgical alveolus and begin the installation by turning the driver clockwise. Complete the implant installation until planned final position. The installation driver must be connected to the AllSelect Prime torque wrench (9780712), or a counter-angle motor. It must be ensured that the torque, during installation of the implant, does not exceed the recommended maximum torque of 70 N.cm.

#### PROSTHETIC SEQUENCE

In the immediate post-surgical period, according to the chosen loading protocol, install the implant cover screw, healing cup or prosthetic component, respecting their indications and limitations.

Use the AllPrime Select height selector (9780768) to set the prosthetic component selection.

#### PRE-OPERATIVE AND PROSTHETIC PLANNING

The model, diameter, length, position and quantity of implants should be selected for each clinical case, taking into consideration the anatomy, region to be rehabilitated, bone quality and quantity, and available space. Diagnostic closure of the clinical case should be performed, ensuring adequate alignment between implant (s), prosthetic component (s) and prosthesis (s). The permitted angulation limit for implant installation is 30°, since it is the maximum angulation present in the AllSelect Prime prosthetic components.

#### PRESENTATION

The AllSelect Prime implants are conditioned in a package that offers triple protection: the primary package that contains the implant and is composed of titanium housing, plastic tube and flexible cap, the secondary package which contains the tube and it's composed of blister sealed with tyvek surgical paper, and the tertiary

package containing the sealed blister and is composed of folded carton paper and sealing by tamper proof label. This external label contains the identification and all information pertinent to the implant to be installed. Inside the tertiary carton, we also have the traceability label that is detachable on 3 smaller labels: one for the user to place on the patient's medical record, another for the patient himself and the last label for reserve.

#### STERILIZATION

This product is indicated for single use and supplied sterile by the gamma radiation method, therefore, they are supplied ready to use.

#### PRECAUTIONS

- Do not use if the product packaging is broken;
- The implant sterilization is only guaranteed if the blister is not damaged;
- This product should be used immediately after opening the blister, at the time of surgery. If you don't use, dispose it;
- Do not use the implant with expired validity;
- Before installing the implant, make sure that it has the appropriate using way for the selected prosthetic system;
- Observe the intraoral tissues conditions, bone quality and bone quantity of the implant surgical bed through clinical, radiographic and / or tomographic exams. Failure to perform pre-surgical evaluation may compromise the success of the procedure;
- About the systemic aspect, consider the general health of the patient. In particular, care should be taken in cases of patients with drug allergies, local or systemic factors that may interfere with the healing processes of bone or soft tissues or in the process of osseointegration, for example: bone already exposed to radiation in the head or neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, treatable maxillary pathologies and oral mucosa alterations;
- Bisphosphonate therapy results in potential risk of peri-implant osseonecrosis;
- The material to be used during the procedure must be sterile;
- Before each procedure, make sure that the parts fit together perfectly;
- The maximum installation torque is 70 N.cm. Insertion torque greater than that maximum, may can make the system inoperative;
- Ensure that the parts are not swallowed or

aspirated by the person who is receiving the product (patient);

- This product is for single use and should not be resterilized;
- The reusing of this product may cause:
  - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
  - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
  - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- When installing components on the implant, be sure to align the connection or the thread, thus avoiding damage inside the implant or locking the system;
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding compromising the implant / prosthesis assembly;
- Before each procedure, check the conditions of the AllSelect Prime instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear;
- Always use the sequence of AllSelect Prime products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllSelect Prime from offering any product warranty;
- Instruct the patient that if he or she notices that the prosthesis will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component;
- The clinical control of each prosthetic work executed with implants should be done periodically, so that the implant dentist can check the balance of the system or intervene clinically, if necessary;
- In the case of prosthetic maintenance where the morse taper component need to be removed, use the AllSelect Prime abutment remover (9780710) for disassembly in order to avoid unnecessary forces to the implant;
- Although the implants or components are

exceptionally resistant, they are subject to fractures in adverse biomechanical situations and / or inadequate rehabilitation plans.

- Do not exceed the recommended torque in the healing cups, implant cover screws, transfer abutments, as they will be removed and may rupture the implant bone interface.

#### **ADVERSE EFFECTS**

The installation of implants, just like any other surgical procedure, can cause mild discomfort and localized edema. More persistent symptoms may occur as: chronic pain related to the dental implant, permanent paraesthesia, dysesthesia, loss of maxillary / mandibular bone border, localized or systemic infection, oro-antral or oro-nasal fistula, adjacent teeth adversely affected, irreversible damage to the adjacent teeth, fractures on implant, maxillary, bone or prosthesis, prosthetic problems, nerve damage, exfoliation, hyperplasia.

The eventual failure of the osseointegration process, or even after osseointegration, with the loss of the implant and the prosthesis can occur due to: inadequate osteotomy, infections, diseases or systemic problems, decrease in remaining bone quantity or quality, parafunction, lack or failure of irrigation, use of non-specific instruments and / or no cutting power, poor oral hygiene, lack of cooperation on the part of the patient, occlusal trauma, lack of prosthetic passivity and absence of specific training.

#### **PÓS-OPERATIVE CARE AND MAINTENANCE**

Orient the patient about the need to perform a professional follow-up after surgery and obey the care guidelines, hygiene of the surgical area and prosthesis, as well as the prescription of medications. These guidelines are the responsibility of the professional.

#### **STORAGE CONDITIONS**

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

#### **DISPOSAL OF MATERIAL**

All products and consumables used in a surgery for the implant dental installation can endanger the health of those who handle them after use. Before disposing them in the environment, it is recommended to consult and comply with the current regulations.

#### **DATE OF EXPIRATION**

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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