## IFU


(Instruction for use)

### Macho de Rosca / Bone Tap / Macho de Roscar

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthetists and Oral Maxillofacial Surgeons.

**DESCRIPTION**

Bone Tap are instrumentally used as an option in dental implant surgery. The operating principle of Thread Screw is based on the formation of threads in the socket, reducing approximately 30% of the installation torque and thus protecting the prosthetic connection of the implant. It has at one end a recess (square 3.9) for connection of the ratchet and at the other end an active portion with cutting edges for thread formation and laser markings regarding the length of the AllPrime Select implants.

All bone tap models are shown in Table 1 below:

Sales Code	Active Diameter	Male Thread Usage Indication		Recommend d $\varnothing$ drill
		Implant		
		CM	HE	
9780703	$\varnothing$ 3,4	$\varnothing$ 3,5	$\varnothing$ 3,5	$\varnothing$ 3,2
9780704	$\varnothing$ 3,9	$\varnothing$ 4,0	$\varnothing$ 4,1	$\varnothing$ 3,8
9780705	$\varnothing$ 4,9	$\varnothing$ 5,0	$\varnothing$ 5,0	$\varnothing$ 4,8

The Bone taps receive the application of a thin biocompatible carbon coating that imparts dark coloration to the cutting surface. This film has as main purposes:

- Increased resistance to oxidation;
- Reduced friction between drill and bone;
- Increased wear resistance;
- Reduction of the heat generated in the bone during the osteotomy;
- Optimization of the visualization of height markings;

The Bone taps are made of tempered stainless steel surgical steel, according to ASTM F899.

**APPLICATION**

The Bone Tap is indicated when the installation torque of the implant exceeds the recommended maximum torque. The threaded male acts on the formation of threads in bone bed type I and II during the surgical procedure, according to the bone classification of Lekholm and Zarb (1985).

Bone classification of Lekholm and Zarb (1985):

- Bone Type I = Homogeneous compact bone.
- Bone Type II = Compact bone with trabecular bone nucleus.
- Bone Type III = Fine layer of cortical bone around dense trabecular bone with good resistance.
- Bone Type IV = Fine layer of cortical bone around trabecular bone with low density.

The AllPrime Select Bone Tap are available in a single length, serving the entire oral region. The length of the AllPrime Select Tap was designed to operate in the same occlusal space as the surgical drill.


**CONTRAINDICATIONS**

Bone Taps should not be used on bone types III and IV.

Do not use Screwdriver with perforations outside the indication (see table 1).

**HANDLING**

For the use of the AllPrime Select Bone Tap, first use the sequence of drills according to the planning and the chosen implant. Use the bi-digital gear, the Torque Meter and capture the Screwdriver with the square (3.9mm) ratchet fitting side. Insert the Screw Into the surgical socket, keeping the drill shaft. Make insertion movements clockwise by lightly pressing the bi-digital gear with the torque wrench. Insert the Screwdriver until the marking for the chosen implant (see image below). There may be a need to change the direction of insertion and removal a few times for thread formation.


The 3.4 mm bone tap should be used after drilling with the  $\varnothing$ 3.2 mm drill, the bone tap 3.9 mm after drilling with the  $\varnothing$ 3.8 mm drill bit and the bone tap 4.9 mm after drilling with the  $\varnothing$ 4.8 mm bit.

To remove the threaded male from the surgical socket, reverse the direction of the bi-digital gear in the torque wrench counterclockwise and proceed with the complete removal of the threaded thread, otherwise the thread shape may be damaged.

**SANITATION AND STERILIZATION**

Screwdrivers are reusable, supplied non-sterile and must be sanitized and sterilized before and after use, following the recommendations of the document available at the link below:

**PRESENTATION**

All AllPrime Select Bone Taps are individually packaged in tyvek® surgical paper packaging and clear film display. The instruments undergo a decontamination process that eliminates any type of dirt, but the hygiene process, described above, must be followed as soon as it is first used.

**PRECAUTIONS**

- Do not use if the packaging has been tampered

with;

- This product should be used sterile. Do not sterilize under dry heat (greenhouse), always use steam autoclave and pressure;
- Chemical sterilization is not recommended;
- Do not autoclave the tapping screw in its original package;
- Avoid lateral loading on the tap;
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient);
- Before each procedure check the conditions of the instruments in fi ls, always respecting their useful life. Replace instruments in case of damage, missing markings, deformation and wear;
- Before each procedure, make sure that the threaded male is properly engaged in the two-wheel gear, otherwise deformation, derailing, or even breaking of the tapping male may occur in the torque application;
- The use of excessive force can compromise the surgical act and damage the male from threading;
- Always use the product sequence In fi ls. The use of instrumentation and / or prosthetic components from other manufacturers does not guarantee the perfect function of the Implant System and exempts any warranty of the product.

**ADVERSE EFFECTS**

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can cause damage to the surgical sequence of the implant to be installed.

**STORAGE CONDITIONS**

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

**DISPOSAL OF MATERIAL**

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

**LIFE CYCLE**

The AllPrime Select hand-driver should be disposed when it loses its functionality.

**SHELF LIFE**

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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