**IFU**

*(Instruction for use)*

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

## DESCRIPTION

The AllPrime Select surgical hammer (9780815) has an anatomical handle for instrument manipulation, and at one end it has a body with an active steel or polymer tip for impact application at the end of the indicated instruments. The connection of the cable with the active body and the body with the polymer tip is through thread, so they can be easily disassembled for hygiene. The raw material used is stainless steel surgical series 300, according to ASTM F899 303, it also has an active tip option produced in polymer.

## APPLICATION

The hammer is a surgical instrument indicated to help in dental procedures. The product is used in surgical procedures for the compacting tips and expander in the application of impact force at the end of the indicated instruments.

## CONTRAINDICATIONS

The surgical hammer does not present contraindications, provided they are used correctly for the purposes indicated.

## HANDLING

Once sterilized, the surgical hammer should be used in compliance with all biosafety and surgical asepsis requirements.

All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the instruments in the best possible way, respecting the AllPrime Select indications for each instrument.

Hold the surgical hammer firmly near the end of the handle manipulation. Whenever hammering, make sure that the active end is parallel to the surface to be hammered. Raise the active tip of the hammer with the metal or polymer surface facing the desired impact surface, and move it quickly and steadily. To disassemble the surgical hammer, firmly hold the handle and rotate the active body counterclockwise until it unscrews the handle. To assemble, run the cable and body joint by screwing the parts clockwise, forcing at the end to ensure that the parts do not come loose. The same process is used to disassemble the polymer tip of the active body.

## SANITATION AND STERILIZATION

The surgical hammer is reusable, provided non-sterile and should be sanitized and sterilized before and after use by following the steps mentioned in the link below:

Surgical Kit Cleaning Manual

## STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

## DISPOSAL OF MATERIAL

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

## LIFE CYCLE

The AllPrime Select surgical hammer should be disposed when it loses its functionality.

## SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

## PRESENTATION

The AllPrime Select surgical hammer is individually packaged in tyvek® surgical paper and transparent film display. The instrument undergoes a decontamination process that eliminates any type of dirtiness, but the sanitation process, described above, must be followed as soon as it's first used.

## PRECAUTIONS

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- Before each procedure, make sure that the connection between the cable and the body is firm;
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be difficult;
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient);
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, thread loosening or fracture of the prosthetic screw;
- Make sure you have all the necessary instruments for the surgery according to surgical planning;
- The use of excessive force can compromise the surgical act;
- Before each procedure, check the conditions of the AllPrime Select instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear;
- Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty.

## ADVERSE EFFECTS

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can cause damage to the surgical sequence of the implant to be installed.

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