IFU




(Instruction for use)

Parafuso / Screw / Tornillo

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The screws are part of the clinical and laboratory family components and have two distinct parts: head and body. The head has a hexagonal internal geometry for connection of the key that will allow the application of the installation torque. The body has thread with geometry that will allow the connection to the specific component.

The raw material used is titanium alloy (grade 5), according to ASTM F136 for clinical screws or 300 series surgical stainless steel, according to ASTM F 899 303 for laboratory screws.

The screws are presented in the table below:

Sales Code	Type		Application	Thread	Max. Torque	Connection
	Clinical	Lab				
9780545	X		Pillars, Abutments, Uclas, Temporary and Angled Standard Base Abutments for MT Implants	M1.6	25 N.cm	Hex. 1,2
9780380		X	Open Tray Transfers and Burnout of Mini 4.0 and 5.0 mm and Angled Standard Base Abutments	M1.4	10 N.cm	Hex. 1,2
9780594 9780595 9780596 9780597	X		Base Abutments Mini Standard 4.8mm for Morse Taper Implants (all heights)	M1.6	25 N.cm	Hex. 1,2
9780379	X		Uclas and Temporary for Base Abutments Mini 4.0 and 5.0 mm and Angled Standard	M1.4	10 N.cm	Hex. 1,2
9780281	X		Pillars, Abutments, Uclas and Temporary for Internal Hex. and External Hex. Implants	M2	30 N.cm	Hex. 1,2
9780605 9780606 9780607	X		Ball Abutment for Morse Taper	M1.6	25 N.cm	Qdr. 2,0
9780334	X		Uclas and Temporary for Multi-Unit and Non-Rotational Base Abutments	M2	30 N.cm	Hex. 1,2
9780332		X	Closed Tray Transfer and Burnout Base Abutments	M2	15 N.cm	Hex. 1,2
9780333		X	Open Tray Transfer and Burnout Base Abutments	M2	15 N.cm	Hex. 1,2
9780381	X		Angled Standard Base Abutment for Internal Hex. and External Hex. Implants	M2	30 N.cm	Hex. 1,2
9780546		X	Open Tray Transfer and Burnout Abutments for Morse Taper Implants	M1.6	10 N.cm	Hex. 1,2
9780282		X	Open Tray Transfer and Burnout Abutments for External Hexagon Implants	M2	15 N.cm	Hex. 1,2
9780547		X	Closed Tray Transfer for Morse Taper Implants	M1.6	10 N.cm	Hex. 1,2
9780283		X	Closed Tray Transfer for External Hexagon Implants	M2	15 N.cm	Hex. 1,2
9780587 9780588 9780589	X		Multi-Unit Base Abutment for Morse Taper Implants (all measures)	M1.6	25 N.cm	Hex. 1,2
9780329 9780330 9780331	X		Multi-Unit Base Abutment for External Hex. Implants (all measures)	M2	30 N.cm	Hex. 1,2

APPLICATION

Clinical and laboratory screws are indicated for definitive or temporary fixation of their respective components according to the recommendation described in the previous table. Before finalizing the clinical work or any manipulation in the laboratory, make sure that the elements have been tightened with the required torque.

CONTRAINDICATIONS

We do not recommend the use of laboratorial screws in the mouth if the patient shows any signs of allergy or hypersensitivity to the chemical components of the material: 300 series surgical stainless steel, according to ASTM F 899 303.

In the case of laboratory screws, there is no

contraindication, provided that its purposes or indications are respected.

MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information | Magnetic Resonance Imaging (MR)” at <https://www.allprimedental.com>

HANDLING

Collect the clinical or laboratory screw using the length of the driver most appropriate to the work you want to perform. At the time of installation, be sure to align the screw along the component’s shaft. The component must be perfectly seated to the implant, transmucosal or analogous after application of the necessary torque on the screw. The maximum torque to be applied to each screw must be respected.

PRESENTATION

Any screw is indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization By Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient);
- Any type of screw is for single use only and should not be re-sterilized.

- Reusing this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- At the time of installation, be sure to align the screw along the component shaft, thus avoiding damage inside the implant or locking the system;
- The component must be perfectly seated to the implant, base abutment or analog after the necessary torque has been applied to the screw. This check is recommended using periapical radiographs;
- It is recommended that the final prosthesis be sanitized before installation in the mouth;
- Instruct the patient that if he or she notices that the healing cup will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the

different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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