IFU

(Instruction for use)

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

AllPrime Select Parallels and depth gauges are supplied in 3 models (9780707, 9780708 and 9780709) and have a center ring with diameters of 3.5mm (9780707), 4.0mm (9780708) and 5.0mm (9780709), which may refer to the distal median space for definition of the most suitable implant or component. This central ring has a small hole to pass and tie a dental cord for safety. At the upper end there is a rod with a diameter of 1.8 mm, the product is marked by laser-engraved grooves that refer to the lengths of implants in the In f is system, so that the user can detect the depth of the drilling performed. The parallelometer and depth gauge is an integral part of the surgical kit In f is and comes in a quantity of 6 pieces (2 pieces of each model) so that the user can compare the parallelism between perforations.

The AllPrime Select parallelometers and depth gauges are made in titanium (grade 5), according to ASTM F136 and anodized for model differentiation.

The model

9780707 is anodized in purple (diameter 3.5), 9780708 in yellow (diameter 4.0) and 9780709 in blue (diameter 5.0). See table 1

Model	Color	Image	Diameter
9780707	Lilac		Ø3,5mm
9780708	Yellow		Ø4,0mm
9780709	Blue		Ø5,0mm

Table 1

APPLICATION

The AllPrime Select parallelometer and depth gauge is a surgical instrument used in dental surgeries to check the angulation and depth by inserting the part into the drilling performed. The parallelometer and depth gauge is available as follows:

The smaller diameter tip of 1.8 mm is used to verify the parallelism after the drilling with the 2.0 mm initial drill, in this verification we can also verify the drilling depth performed.

At the other tip we have a diameter of 2.4 mm, used to verify the parallelism after drilling with the 2.5 mm helical drill and the 2.7 mm diameter, used to verify the parallelism after drilling with the helical drill of 2.8 mm. We use this condition because in the first drillings we can still make small corrections in the parallelism in order to make the positioning of the implants more efficient.

CONTRAINDICATIONS

The parallelometer and depth gauge do not present contraindications, provided they are used correctly for the purposes indicated.

HANDLING

Once sterilized, the parallelometer and depth gauge should be used in compliance with all biosafety and surgical asepsis requirements.


All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the instruments in the best possible way, respecting the AllPrime Select indications for each instrument.

Before using the product, it is recommended to tie a dental floss into the hole to avoid the risk that the patient can swallow or suck the part.

After drilling with the 2.0 mm initial drill, the parallelometer and depth gauge is inserted into the drill by the smallest tip to verify the drill angle and depth. After drilling 2.5 mm and 2.8 mm diameter drills (if applicable), it is recommended to check again the parallelism and, if the drilling has a small angular deviation, we can still correct it with the subsequent drill. This small angular correction is only possible up to a diameter of 3.2 mm.

SANITATION AND STERILIZATION

The parallelometer and depth gauge is reusable, provided non-sterile and should be sanitized and sterilized before and after use by following the steps mentioned in the link below:

PRESENTATION

The parallelometer and depth gauge is individually packaged in tyvek® surgical paper and transparent film display. The instrument undergoes a decontamination process that eliminates any type of dirtiness, but the sanitation process, described above, must be followed as soon as it's first used.

PRECAUTIONS

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be difficult;
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient);
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- Make sure you have all the necessary instruments for the surgery according to surgical
- planning;
- Before each procedure, check the conditions of the AllPrime Select instruments, always respecting their useful life. Replace the instruments in case of damage,

erasure, deformation and excessive wear;

- Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty.

ADVERSE EFFECTS

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can cause damage to the surgical sequence of the implant to be installed.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

LIFE CYCLE

The AllPrime Select hand-driver should be disposed when it loses its functionality.

SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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