IFU

(Instruction for use)

Pilar Reto Cimentado / Cement Straight Abutment / Pilar Reto Cimentado

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The straight abutments for cementation are part of the prosthetic components family and are used as an intermediate between the implant and the prosthesis (crown). The product is available with different connections, depending on the type of implant used.

The abutments have a slightly conical surface on the prosthesis fitting part to facilitate the

excess exit of the cement and ensure a uniform anchorage of the crown.

The raw material used is the titanium alloy (grade 5), according to ASTM F136.

The straight abutments for cementation are presented in the table below:

Sales Code	Connection Type				Implant Diameter	Design		Collar Height	Cimentation Height	Screw
	MT	ZL	EH	IH		Regular	Aesthetic			
9780293			X		3,5		X	---	7,0	M2
9780294			X		4,0		X	---	7,0	M2
9780295			X		5,0		X	---	7,0	M2
9780554	X				3,5-4,0-5,0		X	1,5	4,0	M1.6
9780555	X				3,5-4,0-5,0		X	3,0	4,0	M1.6
9780556	X				3,5-4,0-5,0		X	4,5	4,0	M1.6
9780290			X		3,5		X	---	10,0	M2
9780291			X		4,0		X	---	10,0	M2
9780292			X		5,0		X	---	10,0	M2
9780557	X				3,5-4,0-5,0		X	1,5	6,0	M1.6
9780558	X				3,5-4,0-5,0		X	3,0	6,0	M1.6
9780559	X				3,5-4,0-5,0		X	4,5	6,0	M1.6
9780296			X		3,5	X		---	7,0	M2
9780297			X		4,0	X		---	7,0	M2
9780298			X		5,0	X		---	7,0	M2
9780641		X			3,5	X		---	7,0	M1.6
9780642		X			4,0	X		---	7,0	M1.6
9780643		X			5,0	X		---	7,0	M1.6
9780548	X				3,5-4,0-5,0	X		1,5	4,0	M1.6
9780549	X				3,5-4,0-5,0	X		3,0	4,0	M1.6
9780550	X				3,5-4,0-5,0	X		4,5	4,0	M1.6
9780299			X		3,5	X		---	10,0	M2
9780300			X		4,0	X		---	10,0	M2
9780301			X		5,0	X		---	10,0	M2
9780644		X			3,5	X		---	10,0	M1.6
9780645		X			4,0	X		---	10,0	M1.6
9780646		X			5,0	X		---	10,0	M1.6
9780551	X				3,5-4,0-5,0	X		1,5	6,0	M1.6
9780552	X				3,5-4,0-5,0	X		3,0	6,0	M1.6
9780553	X				3,5-4,0-5,0	X		4,5	6,0	M1.6

APPLICATION

The straight abutments are indicated for cemented unit prostheses on the implant, according to available interocclusal space, existing gingival height and three-dimensional position of the implant. It can be used in immediate or late rehabilitation procedures, in the maxilla or mandible. Before finishing the clinical work, make sure that the screw has been tightened with the necessary torque.

CONTRAINDICATIONS

Contraindicated in cases where the interocclusal space is insufficient and / or the three-dimensional position of the implant is unsatisfactory.

MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information | Magnetic Resonance Imaging (MR)” at <https://www.allprimedental.com>

HANDLING

It can be used in two ways:

- With the appropriate connection, fit the screw inside the abutment, install it over the implant and prepare it after application of the recommended torque on the screw;
- Transfer the implant installed to the working model (plaster). With the appropriate connection, fit the screw inside the selected abutment, position it on the corresponding analogue in the model and prepare it.

The straight abutments should be used in parallel prostheses without angulation. The cementation of the prosthesis is performed in the mouth following the instructions of the cement manufacturer. In this case, it is worth remembering that it is necessary to know the technique to avoid excess cement and chemical burns in the tissues adjacent to the implant or base abutment. Protect screw access during the cementation process, if applicable.

In the mouth, always apply the recommended torque on the screw.

PRESENTATION

The abutments are indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization by Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Ensure that the pieces are not swallowed or aspirated by the person receiving the product (patient);
- Any type of abutment is for single use only and should not be re-sterilized.
- Reusing this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty;
- Special care should be taken during prosthesis design for what the interocclusal space is not insufficient or excessive;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- At the time of installation, be sure to align the screw along the component shaft, thus avoiding damage inside the implant or locking the system;

- The component must be perfectly seated to the implant, base abutment or analog after the necessary torque has been applied to the screw. This check is recommended using periapical radiographs with the parallelism technique, when applicable.
- When using a unitary restoration, it is medically advisable to maintain a minimum foundation height of 4 mm above the gum;
- It is recommended that the final prosthesis be sanitized before installation in the mouth;
- Instruct the patient that if he or she notices that the healing cup will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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