## IFU

(Instruction for use)

Pino de Fixação de Guia Cirúrgica / Surgical Guide Fixation Pin / Pasador de Seguridad de Guia Quirúrgica

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological beings and therefore, to obtain results, it must have been improved to the original standards according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthetists and Oral Maxillofacial Surgeons.

#### DESCRIPTION

The Surgical Guide Fixation Pin is a cylindrical pin, which has a sharp active tip end, has a Stop for use reference on the washer and serves for bone anchoring. The other end has a cylindrical tip with a groove for handling the Surgical Guide Fixing Pin.

The surgical guide fastening pin is made of titanium alloy, according to ASTM F136.

#### APPLICATION

The Guide Guide Surgical Guide pin serves as a fixation bracket for the prototype surgical guide during the implant installation procedure with AllPrime Select Guide technology.

#### CONTRAINDICATIONS

No contraindications, provided it is used correctly for the indicated purposes.

#### HANDLING

Take the Surgical Guide to your mouth and make sure it is seated in the correct position. Insert the 1.3 Guide Drill into the contra-angle. Set up the surgical motor at a speed of 500-800 rpm. Take Drill 1.3 to the washer of the Surgical Guide Clamp Pin. Start the engine and drill with continuous insertion and removal movements with heavy irrigation until the Drill Stop 1.3 touches the washer of the Surgical Guide Clamp Pin. The drilling direction should follow the direction of the washer, without lateral movements. Irrigation can be manual or combined with engine irrigation. After drilling, manually insert the Surgical Guide Fastening Pin until the Stop touches the washer. After all implants are installed, manually remove the Surgical Guide Fastening Pins and Surgical Guide.

NOTE: Do not interrupt engine rotation during drilling.

#### SANITATION AND STERILIZATION

The Surgical Guide Fixation Pin is reusable, supplied non-sterile and should be sanitized and sterilized before and after use, following the recommendations of the document available at the link below:

Novo AllPrime Select's Surgical Kit

#### PRESENTATION

The AllPrime Select guide surgical guide pin is individually packaged in tyvek® surgical paper packaging and clear film display. The instruments undergo a decontamination process that eliminates any type of dirt, but the hygiene process, described above, must be followed as soon as it is first used.

#### PRECAUTIONS

- Do not use if the packaging has been tampered with;
- This product should be used sterile. Do not sterilize under dry heat (greenhouse), always use steam autoclave and pressure;
- Chemical sterilization is not recommended;
- Do not autoclave the explorer probe in its original packaging;
- During the surgical procedure make sure not to exert lateral movements with the surgical instruments, avoiding the loss of the Washer attachment and fracture of the drill.
- Before each procedure check the conditions of the instruments In file, always respecting their useful life. Replace instruments in case of damage, missing markings, deformation and wear;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly resulting in system failure, such as implant loss or fracture, loosening or fracture of prosthetic components and / or screws.
- Always use the product sequence In file. The use of instrumentation and / or prosthetic components from other manufacturers does not guarantee the perfect function of the Implant System and exempts any warranty of the product.

#### ADVERSE EFFECTS

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can lead to unsuccessful perforation and even damage to the patient.

#### STORAGE CONDITIONS

The product should be stored in a dry place in its original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence the original characteristics of the carton and the product.

#### DISPOSAL OF MATERIAL

The operator must always evaluate the condition of the instrumentation after each use and it must be discarded in accordance with current regulations, complying with the environmental requirements and the criterion of potentially dangerous products, and if the instrument is damaged, markings erased, deformed and worn, taking into account the different levels

of contamination, physical risks and selective collection.

#### LIFE CYCLE

Guide Drill Guides should be discarded when they lose their functionality or the user detects excessive wear.

#### SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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