IFU

(Instruction for use)

Saca Implante / Implant Remover / Saca Implante

PORTUGUÊS

ENGLISH

ESPAÑOL

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PRAZO DE VALIDADE

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This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The AllPrime Select implant remover has a single length, capable of attending to all clinical cases and occlusal spaces. At its active tapered end, the product contains a left-hand taper thread, 0.4 mm pitch, which engages the M2.0 and M1.6 internal threaded implant models. At its upper end, the implant sachet has a universal square fit of 3.9 mm, which connects the AllPrime Select solid square wrench (9780751), specially developed for use in implant implants. The raw material used is stainless steel tempered steel, according to ASTM F899.

The AllPrime Select implant remover is shown in the table below:

Sales Code	Initial thread diameter	Active length	Used in threads
9780760	1,4	19	M2,0 / M1,6

APPLICATION

The implant remover is a surgical instrument indicated to be used manually in the removal of poorly positioned or mechanically damaged implants, which are already osseointegrated. It uses the mechanical counter-torque principle (anti-clockwise) to promote the rupture of osseointegration and deinsertion of the implant by the alveolus itself. As the force for implant removal can be high, you should always use the AllPrime Select solid square wrench (9780751). The 9780760 SACA IMPLANTE is indicated for the removal of several connections: external, internal and cone morse, with threads M1,6mm AND M2,0mm.

CONTRAINDICATIONS


It is contraindicated the use of the implant remover in situations of critical bone thickness. The thickness of the bone boards peripheral to the implants that require removal by buccal and lingual loads is one of the limiting factors of the use.

HANDLING


Once sterilized, the implant bag should be used in compliance with all biosafety and surgical asepsis requirements. All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional, who will be responsible for using the instruments in the best possible way, respecting the indications of AllPrime Select for each instrument.

Removal of implants Internal Hexagon and External Hexagon AllPrime Select with thread M2,0:

- Coupling the implant implant on the Square Wrench;


- Install the implant implant (make sure the screw has been removed from the implant);

- Turn the assembly counterclockwise;

- Remove the implant.


Removal of implants Morse Cone Emulsion with thread M1,6:


- Coupling the implant implant on the Square Wrench;

- Wear the internal thread M1,6 with the aid of a Ø2,0mm carbide drill bit;


- Install the implant implant (make sure the screw has been removed from the implant);

- Turn the assembly counterclockwise;

- Remove the implant.

Important points in using the instrument:

- Avoid side loading.
- The solid square wrench has no torque control.
- The initial force for fracture of the bone interface between the osseointegrated implant and the bone is greater than during the procedure.
- Implants remover can accumulate damages during use, it is recommended to use only.

SANITATION

The implant remover is supplied non-sterile and should be sanitized and sterilized before and after use, following the recommendations of the document available at the link below:

[Surgical Kit Cleaning Manual](#)

PRESENTATION

The implant remover is individually packaged in tyvek® surgical paper packaging and transparent film display. The instruments undergo a process of decontamination that eliminates any type of dirt, however the process of hygiene, described above, must be followed. They need sterilization before clinical use.

STERILIZATION

Supplied non-sterile.

Sterilization: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization by Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- A preliminary planning should be performed, evaluating the proximity of the implant to be removed with adjacent dental roots and noble anatomical structures, as this procedure can cause damage to these structures;
- Evaluate the condition of the implant to be removed, as it may be necessary to use trephine instead of the implant remover;
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be

difficulty;

- Ensure that the pieces are not swallowed or aspirated by the person receiving the product (patient);
- Inadequate planning may compromise the performance of the procedure;
- Before each procedure, make sure that the parts are perfectly inserted;
- Make sure you have all the necessary instruments to perform the surgery according to surgical planning;
- This product is for single use only;
- It is recommended to keep the implant remover axial along the axis of the implant to be removed. The inclination of the implant remover during torque application can result in damage and / or fracture of the implant or even the bone;
- It is the responsibility of the professional to use the AllPrime Select products in accordance with the instructions for use.
- The implant remover can accumulate damage during use, SINGLE USE is recommended.
- Reusing this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
 - *Always use the sequence of AllPrime Select products. The use of instruments and / or prosthetic components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty.*

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

LIFE CYCLE

Single use. The product accumulates damage during use.

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