IFU




(Instruction for use)

Sistema O'ring / Ball System / Sistema de Bola

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to AllPrime Select instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The o’ring system consists of the following prosthetic components: a ball abutment, which has different connections, depending on the type of implant used, as well as collars

corresponding to the gingival heights, taking into account the mucosal thickness variations (see table below), mini-female cap (9780493), spacer ring (9780494), spacer capsule (9780495) and burnout mini abutment (9780847), the latter is only used in cases where the abutment requires customization and / or correction of angles. The o’ring system functions as the preparation, installation and retention of the overdenture that is positioned between the ball-type fitting and the mini-female. The miniature female component (9780493), which is properly positioned with the help of the spacer ring and the capsule and cemented in the prosthesis, has an elastic polymer inner ring (9780496) that fits under pressure in the ball abutment and has a retention force necessary to support and retain the prosthesis. Thus, the

prosthesis can be easily removed for sanitizing or any other operation and replaced by the patient himself.

The raw material used in the ball abutment is titanium alloy (grade 5), according to ASTM F136, the mini female is composed of the body, which is also made of titanium alloy (grade 5), according to ASTM F136, and still has an inner elastic polymer ring, a rigid, calcinable polymer is used for the spacer ring and the cap.

The titanium ball abutments of the o’ring system are shown in the table below:

Sales Code	Connection Type			Piece		Implant Diameter	Collar Height	Driver Connection	Screw	Max. Torque
	CM	LZ	HE	Single	Compound					
9780390			X	X		3,5	1,5	Qdr. 2,0	M2	30 N.cm
9780391			X	X		3,5	3,0	Qdr. 2,0	M2	30 N.cm
9780392			X	X		3,5	4,5	Qdr. 2,0	M2	30 N.cm
9780393			X	X		4,0	1,5	Qdr. 2,0	M2	30 N.cm
9780394			X	X		4,0	3,0	Qdr. 2,0	M2	30 N.cm
9780395			X	X		4,0	4,5	Qdr. 2,0	M2	30 N.cm
9780396			X	X		5,0	1,5	Qdr. 2,0	M2	30 N.cm
9780397			X	X		5,0	3,0	Qdr. 2,0	M2	30 N.cm
9780850			X	X		5,0	4,5	Qdr. 2,0	M2	30 N.cm
9780602	X				X	3,5/4,0/5,0	1,5	Qdr. 2,0	M1.6	25 N.cm
9780603	X				X	3,5/4,0/5,0	3,0	Qdr. 2,0	M1.6	25 N.cm
9780604	X				X	3,5/4,0/5,0	4,5	Qdr. 2,0	M1.6	25 N.cm

The other parts that make up the o'ring system, mini-female cap (9780493), spacer ring (9780494), capsule (9780495) and burnout mini abutment (9780847), have no connection variation and connects to any model or implant diameter.

APPLICATION

The o'ring system is indicated to compose the retention assembly in implant-supported overdenture removable prosthesis, planned according to each case, aiming at function and aesthetics. While choosing the o'ring system, make sure that the system is best suited for the job to be done. Before finishing the clinical work, make sure that the screw has been tightened with the necessary torque.

CONTRAINDICATIONS

The o'ring system for overdenture is contraindicated in cases where the angulation between implants or roots exceeds 28°.

MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION

These products are made from a metallic material that can be affected by MR energy.

For more information, see "Safety Information | Magnetic Resonance Imaging (MR)" at <https://www.allprimedental.com>

HANDLING

Titanium Ball Abutments:

Install it over the implant of the same connection, and with the correct height of the collar that meets the corresponding gingival height of the patient. Attach to the ball abutment, the plastic spacer capsule (9780495) and mold it, so that is fixed in the mold. Remove the plastic spacer capsule from the mold and pour up the plaster. Obtained the working model (plaster), conduct the prosthesis preparation according to current techniques.

Once the prosthesis is finished, without the complete retentive system, install the prosthesis on the patient and adjust it, over as many sections as necessary. After adjusting the prosthesis, attach the spacer ring (9780494) on the ball abutment, with its lowered design facing upward, and the mini-female cap (9780493), which will fit into the recess of the ring. Put the prosthesis in position and check if its seating and comfort is adequate. Fill the existing space in the prosthesis (destined to receive the mini-female 9780493) with self-polymerizing resin, in a fluid consistency. Take the prosthesis in position and ask the patient to occlude, keeping the prosthesis in position and awaiting

the resin polymerization. Remove the prosthesis and, if present, remove the excess resin from the mini metal capsule and finish it. Remove the spacer ring from the ball abutment, install the prosthesis, and make your final adjustment.

The spacer ring (9780494) is indicated to maintain the space between the mini-female cap (9780493) and the ball abutment or the implant, compensating for the resilience of the mucosal fiber. It is essential to use the spacer ring in the installation of the o'ring system. Always pay attention to the recommended torque for the ball abutments.

Burnout Mini Abutment:

Indicated to correct divergent and / or convergent implants or roots, creating a parallelism and thus improving the axis of removal and insertion of the prosthesis, increasing the durability of the elastic ring. There is a need to use an abutment with connection to the implant installed for the complete personalization of the burnout mini abutment. The burnout mini o'ring abutment features a rod at the top that serves as an attachment knob in the eyeliner for personalization. This rod must be removed after finding the most appropriate position of the o'ring and before running the micro-casting. The product also has a stem at the bottom for root use, if any, and there is no need for implant installation.

Transfer the installed implants to the working model (plaster).

After personalizing the o'ring mini abutments with the help of an eyeliner, correcting the divergences and / or convergences, respecting the gingival height, and using a abutment to capture the correct connection of the implant so that the prosthesis fits more smoothly and precision, execute the micro-smelting respecting all the needs and requirements that the process demands. With properly molded and cast metal parts, follow the process described above for handling the titanium ball abutments. At the end, use the standard screw to fix the prosthesis corresponding to the model of abutment used in the fitting with the implant, always paying attention to the recommended torque.

Elastic Ring Exchange:

The mini-female component (9780493) has an elastic polymer inner ring (9780496) that snaps into the ball-type socket of the ball abutment and has a retention force required to support and retain the prosthesis. With the frequent removal and insertion of the prosthesis by the patient or even by the lack of parallelism between the fittings, the ring (9780496) can lose its elastic characteristic or even if it ruptures naturally due to material wear. If this

ring needs replacement, remove the prosthesis from the patient's mouth and with the help of a curette, completely remove the elastic ring from the inside of the mini female if it is still there. Sanitize the prosthesis before installing the new elastic ring. Install the new elastic ring by pressing it and positioning it with a curette until it is properly seated in the inner channel of the mini-female.

PRESENTATION

All the parts that make up the o'ring system are indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization By Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 °C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Before installing the implant, make sure that it has the appropriate using way for the selected prosthetic system;
- In situations of severe bone atrophy, it is suggested to carefully evaluate the indication of removable prosthesis, because the stability of the rehabilitation is also related to the support in the mucosal;
- It is recommended to use the spacer ring (9780494) between the ball-type socket and the mini-female cap (9780493) to determine the correct distance between these components;
- The absence of the spacer ring (9780494) can impair the functionality of the prosthesis

and cause premature wear of the elastic ring (9780496), as well as causing collision between the ball and mini-female cap fittings (9780493), leading to breakage of the fitting and disable the restraint system;

- When finishing the prosthesis, care must be taken not to damage the mini-female cap (9780493), if it is already installed in the prosthesis;
- The elastic ring (9780496) inside the mini-female cap (9780493) should be replaced regularly every 8 months or when the function is lost. Orient the patient to appear in the clinic regularly, to check the retention of the prosthesis;
- Ensure that the pieces are not swallowed or aspirated by the person receiving the product (patient);
- Any component of the o'ring system is for single use only.
- Reusing this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Always use the sequence of AllPrime Select products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts AllPrime Select from offering any product warranty;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or o'ring interface, thread loosening or fracture of the prosthetic screw or ball abutment;
- At the time of installation, be sure to align the screw or the ball abutment thread along the component shaft, thus avoiding damage inside the implant or locking the system;
- The ball abutment must be perfectly seated to the implant after the required torque has been applied to the abutment itself or to the screw. This check is recommended using periapical radiographs with the parallelism technique, when applicable.
- It is essential to regular hygiene of the prosthesis by the patient. Inefficient hygiene of this type of removable prosthesis can lead to failure of osseointegration with loss of implant and prosthesis. Therefore, the

professional must instruct and train the patient about the proper maintenance of their prosthesis.

- Instruct the patient that if he or she notices that the healing cup will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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